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
Sample preparation for spectroscopy

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
This series of Lifumat furnaces has been developed to offer many advantages to Metallurgical Research and Quality Control environments within industry when compared with conventional techniques. Whether individual samples or routine controls with up to 250 samples in three shift operation, the variety of

application from platinum to iron scrap, highreprodu-cability, low element losses, high working speed and highest service friendliness is proven due to the multiplicity of satisfied customers worldwide. Their high melting capacity with low energy consumption, easy operation via touchscreen, monitoring of safety devices and a limited space requirement are the most out-standing features of this series. Lifumat-Met-6.6P, - 7.7P, -Super and their va-cuum versions are also suitable for production of calibration standards (cross sectional samples) and scrap (recycling) up to 2 kg.

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


Remelting with a LINN unit by induction heating provides the metallurgical laboratory in research, production and recycling industries with a system for single and routine testing which offers many advantages compared with traditional procedures of sample preparation. The important advantages are **reproducibility**, **low element losses**, as well as **high throughput** which are confirmed by the installation of more than 6000 units for all kinds of application from precious metals to scrap recycling.

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With tilting device. Especially suitable for production of metallic standard samples and for pouring into moulds or ceramic moulds. MF-output power 15 kW, working frequency 10 - 50 kHz or 50 - 100 kHz. Controlled by Simatic S7. Protective gas pressure up to 200 mbar, optionally up to 1.5 bar.

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Melting of jewellery alloys (Gold, Silver, Brass, Copper, Bronze, Steel e.g.) and casting in a mould. Inert gas flushing. Also for granulation. MF-output power 10 kW, nominal working frequencies 10 up to 30 kHz. Controlled by Simatic S7. Lifting device for crucibles. Also tiltable version available.

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

The well established fusion system Lifumat-O-2.0P is available with a touch screen controller and program memory for maximum easy of use. A tilting and turning device (O-2.0P-CRS) improves the homogeneity and can also be used for high melting materials. Fusion temperatures up to 1500 °C will be achieved.

An adjustable cooling system ensures best cooling of the sample.

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High-frequency induction heating

By this method of heating fusion in Pt/Au or graphite crucibles with temperatures up to 1500 °C can be achieved within 3 min. The advantage of **induction fusion** is easy handling and **operation safety** at increased temperature with **reduced time**. The eddy currents guarantee that the melt has an optimal mix by the so-called bath movement. By that a **constant and homogeneous** quality of the melt is achieved. Analysis accuracy and economy with optimal sample quantity confirm the high quality of this method with a minimum of element losses, which can be rectified. Graphite-, glass-carbon-, Zirconia- and Pt/Au-crucibles can be applied. Other fusion methods as e.g. heating by torch or within resistance-heated furnaces show many disadvantages: unequal temperature distribution, poor mixing, environmental problems, high operation costs and far more expenditure of time.

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
Semi-automatic high frequency fusion furnace for spectrometer sample preparation

This induction heated fusion system allows the laboratory to produce quickly and economically up to 35 reproducible fusions of oxidic materials in one step. The melting conditions can be selected individually for each sample. This provides a basis for a constant quality in the following evaluation of the fusion for XRF/RFA, AAS, ICP. The main requirements to fusion are: reproducibility, time saving, costs effective, environmental acceptability, no gas burner, flexibility, safe operation, high temperature, ergonomics. These are fulfilled to a high degree by this system. Special attention was given to a high sample throughput and melting power at low energy consumption, simple operation and supervision through operating unit 6F-Vedo, easy use through touch screen controller and programm memory, low space requirement. The entry of sample preparation parameters and sample identification is easy. All processed samples are documented in a file with achieved temperature if a pyrometer is installed. Also, the crucible materials are not limited to Pt/Au. Graphite, glass-carbon or met. zirconia can also be used. Typical applications are: cement industry, mines, foundries, mineralogical institute etc.

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Induction Fusion and Remelt Furnace

This sample preparation unit enables laboratories to carry out both methods of sample preparation in one unit - **remelting of metals and fusion of oxidic samples for XRF/RFA, AAS, ICP, X-Emission.**

The Lifumat-C-3.3P works with a powerful, safe to operate high frequency generator, coupled with a high frequency power switchover for both methods. Control is done by a touchscreen and programm controller with interface for a master computer.

The sample preparation unit can be switched over between manual and automatic operation. Measurement and temperature control by optical pyrometer.

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
max. weight of casting and samples, sample sizes:

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
Up to 50 user-definable programs for both methods. There-after the working steps will run automatically.

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Laboratory furnaces of series VMK and LM

for pre oxidation of oxidic and half-oxidic samples as well as drying of samples before remelting. Useful volume 1 - 25 l.

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Linn High Therm GmbH Heinrich-Hertz-Platz 1 92275 Eschenfelden